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### 2.2.5. RELATIVE DENSITY

The relative density  $d_{20}^{20}$  of a substance is the ratio of the mass of a certain volume of the substance to the mass of an equal volume of water, both weighed at 20 °C.

Determine the relative density  $d_{20}^{20}$  with precision to the number of decimals prescribed in the monograph, using a pycnometer, a density bottle, a hydrostatic balance or a hydrometer. The thrust of air is disregarded during the weighing; this may introduce an error of 1 unit in the third decimal place.

Two other definitions are commonly used.

The relative density  $d_{20}^{20}$  of a substance is the ratio of the mass of a certain volume of the substance at 20 °C to the mass of an equal volume of water at 4 °C.

The density  $\rho_{20}$  of a substance is the ratio of its mass to its volume at 20 °C. It is expressed in kilograms per cubic metre ( $1 \text{ kg}\cdot\text{m}^{-3} = 10^{-3} \text{ g}\cdot\text{cm}^{-3}$ ).

The numerical relationships between the relative density and the density expressed in kilograms per cubic metre are:

$$\rho_{20} = 998.202 d_{20}^{20} \text{ or } 1.00180 \times 10^{-3} \rho_{20}$$

$$\rho_{20} = 999.972 d_4^{20} \text{ or } 1.00003 \times 10^{-3} \rho_{20}$$

$$\rho_4^{20} = 0.998230 d_{20}^{20}$$

### 2.2.42. DENSITY OF SOLIDS

The density of solids corresponds to their average mass per unit volume and typically is expressed in grams per cubic centimetre ( $\text{g}/\text{cm}^3$ ) although the international unit is the kilogram per cubic meter ( $1 \text{ g}/\text{cm}^3 = 1000 \text{ kg}/\text{m}^3$ ).

Unlike gases and liquids whose density depends only on temperature and pressure, the density of a solid particle also depends on its molecular assembly and therefore varies with the crystal structure and degree of crystallinity.

When a solid particle is amorphous or partially amorphous, its density may further depend upon the history of preparation and treatment.

Therefore, unlike fluids, the densities of two chemically equivalent solids may be different, and this difference reflects a difference in solid-state structure. The density of constituent particles is an important physical characteristic of pharmaceutical powders.

The density of a solid particle can assume different values depending on the method used to measure the volume of the particle. It is useful to distinguish three levels of expression of density:

- **the crystal density** which only includes the solid fraction of the material; the crystal density is also called true density;

- **the particle density** which also includes the volume due to intraparticulate pores,

- **the bulk density** which further includes the interparticulate void volume formed in the powder bed; the bulk density is also called apparent density.

#### I. Crystal density

The crystal density of a substance is the average mass per unit volume, exclusive of all voids that are not a fundamental part of the molecular packing arrangement. It is an intrinsic property of the substance, and hence should be independent of the method of determination. The crystal density can be determined either by calculation or by simple measurement.

A. The calculated crystal density is obtained using crystallographic data (size and composition of the unit cell) of a perfect crystal, from for example X-ray diffraction data, and the molecular mass of the substance.

B. The measured crystal density is the mass to volume ratio after measuring the monocrystal mass and volume.

## II. Particle density

The particle density takes into account both the crystal density and the intraparticulate porosity (sealed and/or open pores). Thus, particle density depends on the value of the volume determined which in turn depends on the method of measurement. The particle density can be determined using one of the two following methods.

A. The *pycnometric density* is determined by measuring the volume occupied by a known mass of powder which is equivalent to the volume of gas displaced by the powder using a gas displacement pycnometer (2.9.23). In pycnometric density measurements, the volume determined includes the volume occupied by open pores; however, it excludes the volume occupied by sealed pores or pores inaccessible to the gas. Due to the high diffusivity of helium, which is the preferred choice of gas, most open pores are accessible to the gas. Therefore, the pycnometric density of a finely milled powder is generally not very different from the crystal density.

B. The *mercury porosimeter density* is also called *granular density*. With this method the volume determined also excludes contributions from sealed pores; however, it includes the volume only from open pores larger than some size limit. This pore size limit or minimal access diameter depends on the maximal mercury intrusion pressure applied during the measurement and under normal operating pressures the mercury does not penetrate the finest pores accessible to helium. Various granular densities can be obtained from one sample since, for each applied mercury intrusion pressure, a density can be determined that corresponds to the pore size limit at that pressure.

## III. Bulk and tapped density

The bulk density of a powder includes the contribution of interparticulate void volume. Hence, the bulk density depends on both the density of powder particles and the space arrangement of particles in the powder bed.

The bulk density of a powder is often very difficult to measure since the slightest disturbance of the bed may result in a new density. Thus, it is essential in reporting bulk density to specify how the determination was made.

A. The *bulk density* is determined by measuring the volume of a known mass of powder, that has been passed through a screen, into a graduated cylinder (2.9.15).

B. The *tapped density* is achieved by mechanically tapping a measuring cylinder containing a powder sample. After observing the initial volume, the cylinder is mechanically tapped, and volume readings are taken until little further volume change is observed (2.9.15).

### 2.9.23. PYCNOMETRIC DENSITY OF SOLIDS

The test for pycnometric density of solids is intended to determine the volume occupied by a known mass of powder by measuring the volume of gas displaced under defined conditions. Hence, its pycnometric density is calculated.

#### Apparatus

The apparatus consists of the following:

- a sealed test cell, with an empty cell volume ( $V_c$ ), connected through a valve to a reference cell, with a reference volume ( $V_r$ ),
- a system capable of pressurising the test cell with the measurement gas until a defined pressure ( $P$ ) indicated by a manometer,
- the system is connected to a source of measurement gas, which is preferably helium, unless another gas is specified <sup>[1]</sup>.

The temperature of the gas pycnometer is between 15 °C and 30 °C and must not vary by more than 2 °C during the course of measurement.

The apparatus is calibrated which means that the volumes ( $V_c$ ) and ( $V_r$ ) are determined, using calibrated, polished steel balls having a total volume (around  $6 \text{ cm}^3$  known to the nearest  $0.001 \text{ cm}^3$ ). The procedure described below is followed in two runs. Firstly, with an empty test cell and secondly with the steel balls placed in the test cell. The volumes ( $V_c$ ) and ( $V_r$ ) are calculated using the equation for the sample volume taking into account that the volume is zero in the first run.

### Method

Weigh the test cell of the pycnometer and record the mass. Fill the test cell with a given mass of powder of the substance to be examined. Seal the test cell in the pycnometer. Remove volatile contaminants in the powder by degassing the powder under a constant purge of gas; occasionally, powders may initially have to be degassed under vacuum. Record the system reference pressure ( $P_r$ ) as indicated by the manometer while the valve that connects the reference cell with the test cell is open. Close the valve to separate the reference cell from the test cell. Pressurise the test cell with the gas to an initial pressure ( $P_i$ ) and record the value obtained. Open the valve to connect the reference cell with the test cell. Record the final pressure ( $P_f$ ). Repeat the measurement sequence for the same powder sample until consecutive measurements of the sample volume ( $V_s$ ) agree to within 0.5 per cent. The sample volume is expressed in  $\text{cm}^3$ . Unload the test cell and measure the final powder mass ( $m$ ) expressed in grams.

### Expression of the results

The sample volume ( $V$ ) is given by the expression:

$$V_s = V_c - \frac{V_r}{\frac{(P_i - P_r)}{(P_f - P_r)} - 1}$$

The density ( $\rho$ ) is given by the equation:

$$\rho = \frac{m}{V_s}$$

[1] If gases other than helium are used, it would not be surprising to obtain values different from those obtained with helium, since the penetration of the gas is dependent on the size of the pore as well as the cross-sectional area of the penetrating molecule. For example, the pycnometric density of porous materials will be overestimated by a measure using nitrogen by comparison with helium.

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