
5. vaja

Analiza nukleotidnih zaporedij in molekulsko kloniranje

5. vaja: Analiza nt zaporedij in mol. kloniranje

Pregled vaje

- **raba kodonov**
 - aplikacija tabel rabe kodonov na heterologne gene
- **restriktijske endonukleaze (restriction endonucleases, restriction enzymes)**
 - iskanje mest (linearna, krožna dsDNA)
- **plazmidne karte (plasmid maps)**
- **verižna reakcija s polimerazo (PCR, polymerase chain reaction)**
 - matrica
 - načrtovanje začetnih oligonukleotidov
 - ohranitev bralnega okvirja

v obliki kviza v
spletni učilnici

Raba kodonov

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Verižna reakcija s polimerazo (PCR)

Potek PCR (poenostavljen primer): vezava začetnih oligonukleotidov

5' -AGGCTC**GCGCGGCTGAGAGACTCGCGTAGAGTCGCGTTAGGCTGAAG**CTGA-3'
3' -TCCGAG**CCGCCGACTCTCTGAGCGCGATCTCAGCGCAATCCGACTTC**GA
CT-5'

matrična dsDNA

 **95°C - denaturacija**
(verigi matrične dsDNA se ločita)

5' -AGGCTC**GCGCGGCTGAGAGACTCGCGCTAGAGTCGCCTAGGCTGAAGCTGA**-3'

3' - TCCGAG**CGCGCCGACTCTCTGAGCGC**ATCTCAGCGCAATCCGACTTCGACT-5'

~55°C (T_a) – vezava začetnih oligonukleotidov
(annealing, prileganje; optimalno T_a izračunamo iz T_m)
[temperaturo hitro znižamo s 95°C, da ne nastane matrična dsDNA]

5' -AGGCTC**GCGCGGCTGAGAGACTCGCGCTAGAGTCGGTTAGGCTGAAG**CTGA-3'
 3' -**CAATCCGACTTC**-5'

smerni oligonukleotid (*forward oligo*)

' -CAATCCGACTTC-5'
protismerni oligonukleotid
(*reverse oligo*)

5' - **GCGCGGCTGAG** - 3'
 3' - TCCGAG**CGCGCCGACTCTTGAGCGCGATCTCAGCGCAATCCGACTTC** GACT - 5'

 (se nadaljuje)

5. vaja: Analiza pt zaporedii in mol. kloniranje

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Potek PCR (poenostavlien primer): polimerizacija

↓ ~72°C - polimerizacija
(ta temperatura je optimalna za delovanje termostabilnih DNA-polimeraz)

5' -AGGCTC**GCGCGGCTGAGAGACTCGCGTAGAGTCGCGTTAGGCTGAAG**CTGA-3'
 CTCAGCG**CAATCCGAC**TTCTTC-5'

smerni oligonukleotid (*forward oligo*)

protismerni oligonukleotid
(*reverse oligo*)

5' - **GCGCGGCTGAGAGACTCG** →
 3' - TCCGAG**CCGGCCGACTCTTGAGCGCGATCTCAGCGCAATCCGACTTC**GACT-5'

1


5' -AGGCTC**GCGCGGCTGAGAGACTCGCGTAGAGTCGCGTTAGGCTGAAG**CTGA-3'
 ..AG**CGCGCCGACTCTCTGAGCGCGATCTCAGCGCAATCCGACTTC**-5'

5' - **GCGCGGCTGAGAGACTCGCGTAGAGTCGCCTAGGCTGAAG**CT.
 3' - TCCGAG**CCGGCCGACTCTCTGAGCGCAGTCAGCGCAATCCGACTTC**GA
 GACT-5'

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Oligo Calc: Oligonucleotide Properties Calculator

Enter Oligonucleotide Sequence Name
OD calculations are for single-stranded DNA or RNA

Sequence Backslash Delimited	Length	Molar Weight	G-C content	Tm	Tm (Basic)	Tm (Adj)	Tm (Nearest neighbor)
<input type="text" value="Backslash Delimited"/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> Rev Complement Strand to 2D tRNA <input type="checkbox"/> 5' modification (any) <input type="checkbox"/> 3' modification (any) <input type="checkbox"/> 10 μM Primer <input type="checkbox"/> 1 mM Salt (Na ⁺)				<input type="checkbox"/> Selected molecule <input type="checkbox"/> tRNA			
<input type="button" value="Calculate"/> <input type="button" value="Swap Strands"/> <input type="button" value="BLAST"/> <input type="button" value="mfold"/>							
Physical Constants Length: <input type="text" value=""/> Molar Weight: <input type="text" value=""/> G-C content: <input type="text" value=""/> Tm: <input type="text" value=""/> °C (Basic) 1 mol of s with an Absorbance of <input type="text" value=""/> at 260 nm 1 mol of s is <input type="text" value=""/> nucleotides & contains <input type="text" value=""/> micrograms. Tm: <input type="text" value=""/> °C (Adj Adjusted) Tm: <input type="text" value=""/> °C (Nearest neighbor)							
Thermodynamic Constants/Conditions T = <input type="text" value=""/> K; NaCl at 20°C; pH <input type="text" value=""/> ; μ = <input type="text" value=""/> M							
Rev/RK	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Wt/Wt	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Computed Statistical Properties Calculations 1. <input type="checkbox"/> Minimum base pair required for single primer self-dimerization 2. <input type="checkbox"/> Minimum base pairs required for primer hairpins							
<input type="checkbox"/> Check Self-Complementary							

[http://
www.basic.northwestern.edu
u/biotools/oligocalc.html](http://www.basic.northwestern.edu/biotools/oligocalc.html)

Primer3web version 4.0.0 - Pick primers from a DNA sequence.

disclaimer	code of conduct
Select the Task for primer selection: <input type="text" value="generic"/> <input style="width: 20px; height: 15px; vertical-align: middle;" type="button" value="..."/>	
Paste source sequence below (3'-5'': string of ACCTTNaacgt -- other letters treated as N -- numbers and blanks ignored). FASTA format ok. Please N-cut undesirable sequence (vector, ALUs, LINEs, etc.) or use a Mispriming library (repeat library) : <input type="text" value="NONE"/>	
<input checked="" type="checkbox"/> Pick left primer, or use left primer below <input type="checkbox"/> Pick hybridization probe (internal oligo), or use oligo below <input checked="" type="checkbox"/> Pick right primer, or use right primer (3' to 5' on opposite strand)	
<input style="width: 150px; height: 30px; vertical-align: middle;" type="button" value="Get Primers"/>	

<http://primer3.ut.ee>

Primer-BLAST A tool for finding specific primers.

► NCBI Primer-BLAST: Finding primers specific to your PCR template (using Primer3 and BLAST).

[Reset page](#) [Save search parameters](#) [Retrieve results](#)

PCR Template

Enter accession, gi, or FASTA sequence (A reference record is preferred) [Clear](#)

Or, upload FASTA file [Choose File](#) | no file selected

Primer Parameters

Use my own forward primer
(5'-3' > plus strand)
Use my own reverse primer
(3'-5' > minus strand)

PCR product size
of primers to return

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